

# Spinecare Mom Help Center and Policy Documents

Generated: 2026-06-17

This file is for project archive. Medical, legal, and local-language review is recommended before release.

## User Manual

### Quick Start

After login, complete the profile wizard: nickname, gender, birth date, Cobb angle, curve location, Risser sign, prescribed wearing hours, brace type, and first visit date. The profile is synced to the cloud for the home page, reports, and AI interpretation.

### Bluetooth Device Binding

Open Settings, enter Bluetooth Connection, and scan for devices whose names start with DR-Z-T0. After a successful connection, the app saves the device and will try to reconnect to it next time.

### Read Wearing Data

After Bluetooth connects, the app automatically reads wearing data from the device and uploads it to the cloud. It shows Reading device data, please wait, and then Data reading complete.

### Home Data Insights

The 30-day wearing summary and 7-day wearing chart on the home page use real cloud data. The target standard comes from the prescribed wearing hours in the profile, and colors indicate risk levels.

### Record Features

Skin records are used for photos and notes when a problem is found. Growth records store height. Imaging archives save photos taken by the camera or selected from the phone gallery. These records are included in visit reports.

### Reports and Archive

AI reports, visit reports, and archives are generated from cloud wearing data, skin records, growth records, imaging records, and profile information. An archive keeps a snapshot from the generation time.

### Export and Delete

Export data before deletion as a backup. Before deleting all data, export a backup, tick the confirmation items, enter the child nickname and confirmation text, and wait through the cooling-off period.

## FAQ

### Why bind a Bluetooth device first?

Wearing data stored in the brace sensor must be read through Bluetooth. After binding, the app can connect automatically and upload wearing records.

### What if Bluetooth cannot connect?

Make sure phone Bluetooth is on, the device is nearby, and the name starts with DR-Z-T0. If the last device is not found within 20 seconds, scan and connect again in Settings.

### How is the home compliance rate calculated?

Compliance rate is the percentage of days in the selected period that meet the prescribed wearing hours. The prescribed hours come from the profile.

### Can deleted data be restored?

No. After deletion, cloud health data cannot be restored. To avoid mistakes, the app uses backup check, confirmation text, long press, and a 24-hour cooling-off period.

### Can AI advice replace doctors?

No. AI interpretation only organizes data and provides health education. It does not replace diagnosis, treatment, or brace adjustment by professionals.

# User Agreement

## 1. Service Description

Spinecare Mom helps guardians record and view brace wearing, skin, growth, imaging, and report information. The app provides data organization, reminders, and assisted interpretation.

## 2. Guardian Responsibility

The app is intended for minor health management scenarios. Guardians should confirm they are authorized to enter and manage relevant data and keep the information true, accurate, and current.

## 3. Medical Disclaimer

The app is not a medical diagnosis tool. AI interpretation, report summaries, and reminders are for health management reference only and do not replace medical diagnosis, prescription, follow-up plans, or brace specialist advice.

## 4. Data and Reports

Uploaded or synchronized data may be used for home statistics, intelligent interpretation, visit reports, archive, and export backup. Archived reports keep a snapshot at the time of generation.

## 5. Termination and Deletion

Users may end use through data export and Delete All Data. After deletion, cloud data cannot be restored; please keep a backup in advance.

# Privacy Policy

## 1. Data Collected

The app may collect profile information, wearing records, Bluetooth device information, skin photos and notes, growth records, imaging files, AI consultation content, report archives, and alert messages.

## 2. Purposes

Data is used to show wearing trends, judge compliance, generate visit reports, provide intelligent interpretation, keep archives, export backups, and perform deletion workflows.

## 3. Permissions

Bluetooth permissions connect brace devices. Camera and gallery permissions record skin or imaging materials. Network permission synchronizes data with the cloud database.

## 4. Minor Protection

Minor data should be managed by guardians. The app minimizes unnecessary identity collection and does not require real names or ID numbers by default.

## 5. Delete All Data

Deletion removes the current child's cloud profile, wearing data, records, reports, alerts, device binding, and uploaded files. It cannot be restored after completion.

# Contact Us



## WeChat Support

WeChat Support: Spinecare Mom. Scan the QR code above to add WeChat support.

## Email

Email: [zclei@vip.sina.com](mailto:zclei@vip.sina.com)

## Support Tips

When contacting support, include the child nickname, device name, issue time, and related screenshots or photos. For urgent health issues, contact a doctor or brace specialist promptly.