

Spinecare Mom FAQ

Generated: 2026-06-17

This file is for project archive. Medical, legal, and local-language review is recommended before release.

FAQ

Why bind a Bluetooth device first?

Wearing data stored in the brace sensor must be read through Bluetooth. After binding, the app can connect automatically and upload wearing records.

What if Bluetooth cannot connect?

Make sure phone Bluetooth is on, the device is nearby, and the name starts with DR-Z-T0. If the last device is not found within 20 seconds, scan and connect again in Settings.

How is the home compliance rate calculated?

Compliance rate is the percentage of days in the selected period that meet the prescribed wearing hours. The prescribed hours come from the profile.

Can deleted data be restored?

No. After deletion, cloud health data cannot be restored. To avoid mistakes, the app uses backup check, confirmation text, long press, and a 24-hour cooling-off period.

Can AI advice replace doctors?

No. AI interpretation only organizes data and provides health education. It does not replace diagnosis, treatment, or brace adjustment by professionals.