

Spinecare Mom User Manual

Generated: 2026-06-17

This file is for project archive. Medical, legal, and local-language review is recommended before release.

User Manual

Quick Start

After login, complete the profile wizard: nickname, gender, birth date, Cobb angle, curve location, Risser sign, prescribed wearing hours, brace type, and first visit date. The profile is synced to the cloud for the home page, reports, and AI interpretation.

Bluetooth Device Binding

Open Settings, enter Bluetooth Connection, and scan for devices whose names start with DR-Z-T0. After a successful connection, the app saves the device and will try to reconnect to it next time.

Read Wearing Data

After Bluetooth connects, the app automatically reads wearing data from the device and uploads it to the cloud. It shows Reading device data, please wait, and then Data reading complete.

Home Data Insights

The 30-day wearing summary and 7-day wearing chart on the home page use real cloud data. The target standard comes from the prescribed wearing hours in the profile, and colors indicate risk levels.

Record Features

Skin records are used for photos and notes when a problem is found. Growth records store height. Imaging archives save photos taken by the camera or selected from the phone gallery. These records are included in visit reports.

Reports and Archive

AI reports, visit reports, and archives are generated from cloud wearing data, skin records, growth records, imaging records, and profile information. An archive keeps a snapshot from the generation time.

Export and Delete

Export data before deletion as a backup. Before deleting all data, export a backup, tick the confirmation items, enter the child nickname and confirmation text, and wait through the cooling-off period.